

Uni Bloc

www.shimadzu.com

AUW220D

Max 220g

Min 0.1g

e=0.001g

d1=0.01mg(0-82g)

d2=0.1mg(0-220g)

0.30722 g

SHIMADZU

POWER

BRK

CAL

MENU

→C/T←



UNIT



PRINT





Feb 14, 2026, 16:06



Feb 14, 2026, 16:23

14. Feb. 2026.

Saturday

MoS₂ preparation

Date / /

MoO₃ (5mg)

S (300mg)

0 - 550

110 min

0 - 100

20 min

(80 min Hold)

550 - 650

20 min.

100 - 200

20 min.

650 - 750

20 min.

30 min Hold.

↓

30 (144 min)

↓ cool

34 min

Sublime ⇒ 600 - 700.

~ 160 °C

This time we kept 10 min discrepancy
(10 min shorter to S) ⇒ S will sublime
10 min faster, bcz to make MoS₂,
Sulphur has to reach the substrate
first.

[Regulator repair took a lot of time.]